

#### Use

Corn oil is used for cooking, also as a salad oil and in the manufacturing of margarine.

#### **METHOD 1**

## **Test Equipment**

Instrument: Viscometer or Rheometer

Torque Range: LV

Spindle: YULA-15E

Accessories: Enhanced UL Adapter,

TC-550AP Programmable Bath

Speed: 2, 4, 6, 8 and 10 RPM



#### **Test Method**

The Enhanced UL Adapter<sup>™</sup> may be used with various Brookfield Viscometers or Rheometers. In our example, we used a Brookfield DV3TLV Rheometer with Rheocalc<sup>™</sup> software for automated instrument control and data acquisition. Our test temperature of 25°C was maintained by direct immersion of the Enhanced UL Adapter sample chamber into a TC-550AP Programmable Refrigerated Bath. Representative data from the analyses in Figure 1 show a relatively constant viscosity at 25°C (55 cP to 52 cP) over the speed range tested.




Figure 1: Viscosity vs. Rotational Speed

SETTING THE STANDARDS in Viscosity Testing

#### **METHOD 2**

# **Test Equipment**

Instrument: Viscometer or Rheometer

Torque Range: LV Spindle: LV-1

Accessories: TC-550AP Programmable Bath

Speed: 25, 50, 75 and 100 RPM

### **Test Method**

We used a Brookfield DV3TLV Rheometer with LV guard leg and Rheocalc<sup>™</sup> software for automated instrument control and data acquisition. The sample, tested in a 600 mL low form Griffin Beaker, was immersed in a TC-502P Programmable Refrigerated Bath and equilibrated at 25°C. Representative data from the analyses in Figure 2 show a relatively constant viscosity at 25°C of approximately 52 cP over the speed range tested.




Figure 2: Viscosity vs. Rotational Speed